

P.E.T.S Microchip & Registry System





## Willagee Central LPO PO Box 1033 WA 6156

tel: 1300 537 140 fax: 1300 537 141 acn: 138 062 316

e: sales@homesafeid.com w: www.homesafeid.com

## **ZS007** Animal ID Microchip



- ZS007 is a kind of syringe that designed exclusively for animals. It contains a disposable syringe, a parylene coating microchip, 6 one-dimensional bar code labels and a sterilized paper-plastic packing pouch.
- This product meets all international standards of animal identification (ISO 11784/ISO 11785). Each syringe is produced using top quality European chips with read/write capability and features unique data security. Coupled with Fofia Tech reader series, ZSO07 provides a beneficial solution in animal breeding, slaughter management and animal identification.
- As we adopt the 1.4\*8mm microchip electronic tag produced by the world's leading technology, this syringe is light, handy and more humane to use

Microchip Specification	
R/W standard	ISO11784/5 FDX-B
frequency	134.2KHz
dimension	Ø1.4*8mm
anti-collision standard	IP67
Guarantee	10years
storage temperature	-20°C to +50°C
operation temperature	0°C to 50°C
	Transponder Syringe Specification
color	white
weight	6.86g
length	111±2mm
width	50±1mm
Others	
injection material	PP
packing manner	20 pcs/bag
sterilization	EO gas, 5 years' guarantee
Bar code label	6 adhesive-stickers
International certification	International Committee For Animal Recording (ICAR)

Your pet cannot tell anyone where it lives...HomeSafeID can

HomeSafeID (HSID) unique microchips have a proven record for quality and reliability for use in the electronic identification of companion animals, horses, birds, small animals and wildlife.

HSID microchips are manufactured under the strictest quality control of ISO 9000. All components used in the HSID microchip must be inspected under the AQL 0.65%. Syringe parts must be inspected to a level of 100% viability before releasing to the production line.

## **HomeSafeID Difference**

HSID microchips are delivered direct from the factory to the customer, therefore no distributor costs are imposed which means lower cost for our customers and more timely delivery time frames.

Upon delivery, every single HSID microchip is traceable back to the organisation which it was sold <u>prior</u> to registration taking place.

Because of this unique HSID tracking system transcription errors and duplicates and very rare when used in conjunction with the HSID National Registry. Customers can track chip stock by referring to the CSV file provided with each delivery

This tracking system is a vital safety-net and part of the complete and unique service HSID offers.

## Australian Compliance:

HSID chips conform to AS 5018 which mirrors the International Standard 11784.

The HSID registry complies with all companion animal regulations as mandated by various Australian States and Territories.