

HOMESAFE

P.E.T.S Microchip & Registry System



Willagee Central LPO PO Box 1033 WA 6156

tel: 1300 537 140 fax: 1300 537 141 acn: 138 062 316

sales@homesafeid.com www.homesafeid.com



 ZS008 Animal ID microchip is exclusively designed for animal identification. lt contains disposable needle package, reusable injection gun, a parylene coating microchip, 3 bar code labels and a packing bag pouch. This product is sterilized by EO gas with 5 years guaranteed.

ZS008 Animal ID Microchip

Microchip Specification	
R/W standard	ISO11784/5 FDX-B
frequency	134.2KHz
Dimension	Ø2.12*12mm
anti-collision standard	IP67
Guarantee	10years
storage temperature	-20°C to +50°C
operation temperature	0°C to 50°C
Transponder applicator Specification	
color	blue
weight	46±1g
length	109*76mm
thickness	19mm
	Others
injection material	PP
packing manner	Needle package: 20 pcs / bag Transponder applicator: 1 pc/box
sterilization	EO gas, 5 years' guarantee
bar code label	3 adhesive-stickers
International certification	International Committee For Animal Recording (ICAR)

companion animal regulations as mandated by various Australian States and Your pet cannot tell anyone where it lives...HomeSafeID can

HomeSafeID (HSID) unique microchips have a proven record for quality and reliability for use in the electronic identification of companion animals, horses, birds, small animals and wildlife.

HSID microchips manufactured under the strictest quality control of ISO 9000. All components used in the HSID microchip must be inspected under the AQL 0.65%. Syringe parts must be inspected to a level of 100% viability before releasing to the production line.

HomeSafeID Difference

HSID microchips are delivered direct from the factory to the customer, therefore no distributor costs are imposed which means lower cost for our customers and more timely delivery time frames.

Upon delivery, every single HSID microchip is traceable back to the organisation which it was sold prior to registration taking place.

Because of this unique HSID tracking system transcription errors and duplicates and very rare when used in conjunction with the HSID National Registry. Customers can track chip stock by referring to the CSV file provided with each delivery

This tracking system is a vital safety-net and part of the complete and unique service **HSID** offers.

Australian Compliance:

HSID chips conform to AS 5018 which mirrors the International Standard 11784.

The HSID registry complies with Territories.